

LAB 1 - USE-CASE FLOW

Remote Patient Vitals Alert & Monitoring App

Use Case ID	UC-01
Use Case Name	Record and Monitor Vital Signs
Primary Actor	Remote Patient
Supporting Actor	Vital Signs Device
Secondary Actor	On-Call Caregiver

Main Success Scenario

1. Remote Patient opens the monitoring application.
2. System displays the patient's monitoring dashboard.
3. Patient records the current vital readings, or the Vital Signs Device sends the readings.
4. System validates the received readings.
5. System stores the readings with the patient ID and timestamp.
6. System compares the readings against configured thresholds.
7. System determines that the readings are within acceptable limits.
8. System displays the latest readings and their normal status.
9. Use case ends successfully.

Alternate Flow A - Invalid Vital Reading

- 4a1. System detects that one or more vital readings are invalid.
- 4a2. System displays an error message identifying the invalid reading.
- 4a3. Patient or device is prompted to submit a valid reading.
- 4a4. The use case resumes at Step 3.

Alternate Flow B - Abnormal Vital Reading

- 7a1. System identifies that a vital reading exceeds its configured threshold.
- 7a2. System marks the reading as abnormal.
- 7a3. System generates an alert.
- 7a4. System sends the alert to the On-Call Caregiver.
- 7a5. The caregiver views the alert and may acknowledge it.